

Fig 1: 3D fusion image of cerebral perfusion and angiography (4.6 mSv) (Case courtesy Charité Berlin, Germany)

AquilionTMONE is the first CT scanner capable of imaging whole organ regions up to a width of 16 cm in one rotation and within a split second. Based on the raw volume data, rapid dynamic processes within an entire organ (e.g. heart, pancreas, kidney or brain) may be diagnosed with a time interval of 50 ms, i.e. with a rate of 20 volumes per second. With their smallest effective width of 0.5 mm the detector elements ensure best possible spatial resolution. Image reconstruction of 2 x 320 slices, coupled with special mathematical interpolation, provides a geometric resolution of 0.4 mm image voxels in all directions. Unlike in flying-focus spot

CT: Optimising dosage

By **Dr Jörg Blobel**, Ph.D., Chief Clinical Science, CT Systems Division, Toshiba Medical Systems Corporation, and **Jürgen Mews**, of the CT Systems Division at Toshiba Medical Systems in Neuss, Germany



Dr Jörg Blobel

technology, the signal strengths are not halved and spreadover both slice series. Consequently, this halved signal strength per slice does not have to be compensated by increasing the exposure.

Heart rate adapted temporal resolution between 50 ms and 175 ms allows scanning of the entire organ within one heart beat. For rather high heart rates the temporal resolution of 50 ms is almost half that of a

and non-calcified plaques, the latter frequently resulting in myocardial infarction, can be visualized down to a diameter of less than 1 mm.

This outstanding low-contrast resolution of the Aquilion ONE allows low radiation energy levels at 80 kV, e.g. for diagnostic organ perfusion scans. For the first time perfusion of the entire cerebral volume can be captured simultaneously. The 15-20 volume data sets generated during one minute of the study are acquired with a regular effective patient dose of just 4-6

mSv. Volumetric acquisition permits anatomically accurate fusion of the CT angiography and perfusion volumes (Fig1). The innovative range of CT studies opens up new perspectives for functional diagnostic workup, e.g. in joint movement, peristalsis, dynamic blood flow analysis and perfusion of numerous organs.


In diagnostic pediatric scans the radiation field and the number of detector element rows is collimated to the size of the organ. This rapid examination avoids the need for

breathhold in infants and small children otherwise required by the longer scanning times. If an emergency thoracic CT study is needed in an infant the effective dose at a tube voltage of 80 kV is minimal at 0.16 mSv. Large areas of the body in combination with cardiac CT can be captured by individual volume scans much more rapidly than with helical CT and are subsequently "stitched" to one patient volume. Complex studies of the heart, lungs and the head may be combined with supplementary CT angiography and – based on

the study plan – are optimized by automatic system selection of the scan parameters to reduce radiation exposure. Before the start of the study the expected patient radiation exposure is displayed and can be controlled.

Conventional helical CT mode with the selection of 16, 32 and 64 detector rows is also available and the volume scan mode is enhanced by numerous new dynamic function options. These innovations of the Aquilion ONE will alter the patient workflow between the diagnostic imaging modalities in radiology.

dual source CT scanner and the heart is better locked into position during motion. Compared with a Helical CT unit where the heart volume is captured by individual overlapping rotations the volume scan reduces the effective dose of the normal patient to 1.5-6 mSv – a reduction of 60-80%. Also the disadvantages of the step-and-shoot mode of multislice CT with a scanned field width of 20-40 mm, i.e. the common stepping artifacts at the volume borders and dose doubling along the volume overlap, are overcome. Data capture with Aquilion ONE requires just one heart beat. The modality is robust and offers increased potential to study arrhythmia patients. Calcified



Dose 'hysteria' and the 320-slice CT scanner

For the last few months, **Dr Patrik Rogalla**, Senior Consultant at the Charité, Berlin, specialist in diagnostic radiology and Head of the Computed Tomography Department at Campus Mitte, University Medicine Berlin, has been using a 320-slice CT scanner, one of four currently manufactured. During a *European Hospital* interview, we asked whether the system has provided greater detection rates, and why it has also ignited further debate on radiation dose



marginal. The parameter setting of the equipment, determined by the radiologist based on the individual requirements of a patient, and in particular the clinical question at hand, is much more important.

'It is a little regrettable that the discussion about CT tends to be reduced to a discussion around dose, although the use-risk ratio of CT scanning across medical indications is calculated at around

200:1. This is not only damaging for this technology, but also for radiology itself. This discussion then questions many radiological procedures, such as conventional angiography and fluoroscopy. How can you convince a female patient to have embolisation of a fibroid carried out using fluoroscopy when there is this constant debate over the dangers of radiation during CT scanning? We are taking away our own basis for essential methods of examination and treatment in radiology. I think this is short-sighted. The discussion around dose is necessary, but should be carried out based on the highest levels of knowledge and seriousness. One of the great difficulties with this discussion is that in today's scientific world,

Sometimes radiation dose is a topic discussed hysterically. Dose discussion is necessary, but should be carried out based on the highest levels of knowledge and seriousness

where you cannot refer back to publications from Medline going back more than 4-5 years, nobody can be bothered to carry out proper research into data on exposure to radiation in the way it is required. A piece of research carried out two years ago is referred to because it appeared in a "respectable" journal, and is used as proof that radiation has a "calculable" risk. However, when you take a closer look at this piece of research you then see that it, in turn, refers to another piece of research carried out a few years previously, and this in turn quotes another piece of older research, which eventually takes us back to the original source – data from those who survived Hiroshima. Apart from very few exceptions, these survivors' data are thus the only source – and of course they don't constitute a scientific study

but consist of estimates and observations – observations of a disastrous situation following the dropping of an atomic bomb. This source – and luckily we have not had to experience a more recent, comparable event – is used to deduct a seemingly high degree of evidence, recalculated over and over again and used as hard data to calculate literally the deaths per examination rate, although the basic data hasn't changed and the linearity of radiation damage to dose remains a hypothesis. No new drug, no new examination procedure would ever become market-ready these days if the risk evaluation was carried out without such a lack of solid, scientific base as the discussion around dose in radiology.'

Dr Rogalla is currently using a *dynamic volume* 320-slice scanner. 'This refers to the term dynamic imaging in general,' he explained. 'Dynamics is a process that stretches over a period of time and as such can be captured. This can be movement, perfusion, blood flow or a metabolic process. Previously CT diagnostics was more or less like picture taken with flashlight; the freezing of a moment. However, dynamic imaging captures the course of a process.'

'If you handle dynamic CT imaging wrongly then dose can really become a problem. But we have enough ways around this to avoid that problem. You can distribute the dose across a period of time, literally fraction it. We line up several images one after the other, each of which obviously only requires a fraction of the whole dose, so that the entire dose during an examination of several seconds or minutes remains the same. We compensate for the resulting image loss with mathematical solutions – by averaging over time using new filters. This works. It involves one of the concepts of post-processing; admittedly, this is still under development. I am convinced that we will see very interesting new approaches to solutions during the next scientific congresses.'

'Radiation dose is a topic that is sometimes discussed rather hysterically, and often by many without sufficient basic knowledge on the subject. Reason number one: if people are short on interesting and innovative topics they will talk about radiation dose. This always goes down well and scores points, because radiation is generally perceived as something very negative in our society. Once you tackle this topic you are guaranteed to remain a talking point. The second reason: Of course we must cut down on dose and are obliged by the legislators to examine using the lowest possible dose. The reason why the 320-slice is re-igniting discussions is the potentially lower dose requirement of this new type of CT scanner, and this puts us into the focus of the discussions around dose.'

'There are not quite so trivial discussions around the question of how much dose is required for what kind of image quality. The dose can be measured – it is not that easy but technically it can be done in a fairly reliable way. However, image quality can only be conditionally objectively measured, because a large part of image quality is determined by subjective impression. The real measure should be: How much of a dose do I need for a sensible image quality that allows safe diagnosis? However, seeing how images tend to make very subjective impressions, this opens the floodgates to all kinds of speculation and marketing. CT scanners manufactured by all the main suppliers all require a certain degree of a radiation dose for comparable examinations. The differences between the different manufacturers are rather

MICRODOSE MAMMOGRAPHY. NO IMAGE DISTORTIONS TO HAUNT YOU

Ghosts belong in creaky old buildings. Not in X-rays. Our patented photon-counting technology effectively banishes them. The detector is fast enough to be ready when the next photon arrives, so there's nothing to interfere with image interpretation. For more on the highest image quality at half the radiation dose, visit sectra.com/medical

RIS/PACS | MAMMOGRAPHY | ORTHOPEDICS

SECTRA